

Features

- Compatible with the Mondial® Transporter
 - Converts to chair or stretcher
 - Anodized tubing prevents oxidation for clean hands and clothing
 - Infinitely-adjustable gas-assist backrest, adjusts up to 90°
 - Two wheels (100 mm / 4 in.)
 - Two telescoping handles
 - Two folding handles (foot of stretcher)
 - Folding handles (Chair option)
 - Reversible patient loading: head first or feet first
 - Pivoting footrest is standard
 - Choice of external or internal bolster mattress (sold separately)
 - Optional extendable headrest
 - Swing-down sidearms
 - Features contour shock position
 - Backrest Airway Management Frame (Optional)
 - Includes one 4-point metal buckle restraint and two 180 cm metal buckle restraints
- 60-4300-001 CHAIR-STRETCHER RS-5 (PACKED)
 - 60-4300-001/AW CHAIR-STRETCHER RS-5 WITH FOLDING HEADREST (PACKED)

Specifications

Length	Width	Bed Height	Chair Height	Seat Height	Chair Depth	Seat Depth	Weight	Load Limit
1940 mm	555 mm	160 mm	1390 mm	500 mm	720 mm	410 mm	23 kg	250 kg

Accessories

Headrest Extension

This accessory allows the backrest to be extended and locked into positions up to 18 cm accommodating even the tallest patient

65-4338-001



Backrest Airway Management Frame

The head section of the backrest can be adjusted in three locking positions allowing access to the patient's airway

65-4332-001



Mattress

This 50 mm black vinyl mattress is divided into five sections to contour to our stretcher providing the patient with comfort while being water and stain proof

65-0358-005



IV Pole


This \varnothing h height adjustment which allows operators to set their desired height for treatment

65-1500-030



Set of three patient restraints: one harness torso restraint and two, one-piece belt restraints

63-4400-001



PackRack # 60-1100-001



RS5 Stretcher has been designed to comply with and exceed international standards such as EN 1865 and EN 1789.

